

Plexr[®] OMNIA




***Be the artist
of beauty***

G|M|V
MEDICAL INNOVATION

What

is Plexr®

Omnia 
skin System?

Main indications:

Sagging and excess skin

Subcutaneous fat zones

General signs of skin aging

Points

● ● ● ● ● ● of strength:

Treatment without surgery and stitches.

The nature of the treatments as well as the results are extremely natural with no external substances or fillers are injected.

The Omnia procedures provide a very effective aesthetic solution with full respect for safety and long term health of the skin.

Patient keeps his/her own facial expression, but with a clearly rejuvenated and naturally lifted skin.

The complete Omnia system consists of 3 techniques:

*Plasma Exeresis –
Skin reduction with Plexr®*

*Collagen Shaping with
Needle Shaping*

Lipo Plexr with OFF

*Depending on the chosen treatment,
these are generally synergistically combined but can be applied individually as well.*

Created by nature



Designed by you

Shaped with Plexr Omnia

How to **perform** treatments?

Plexr Omnia Skin System can be performed with Plexr Vecron - Four different technologies in just one medical device. Plexr Vecron includes:

Plexr®

Needle Shaping

OFF

Carboxytherapy



Perform Plexr Omnia

Optional delivery

Do you already have Plexr®? Choose Vecron: a modular medical device – Up to 3 different technologies in the same device: Needle Shaping, OFF, Carboxytherapy



Plexr®

skin reduction

Plexr® is the only patented plasma technology for Aesthetic Medicine.

This specific plasma and its different proprietary impulses have a 3 fold effect on the skin:

Sublimation of the treated skin, resulting in optimal shrinkage and reduction of the skin surface.

In addition, a gentle but energetic bio-boost is given to the deep layers of the skin which results in contraction of the elastic fibres, causing an immediate lifting and stimulation of the skin.

The bio-stimulative impulse to the deeper skin layers ensures an intense and long-lasting build-up of new collagen of type 1, the so-called baby collagen that provides elasticity and a young appearance of the skin. An anesthetic cream is applied locally during this treatment.

Indications for Plexr® treatments among others include: sagging and excess skin, upper and lower eyelid correction, crow's feet (upper lip), wrinkles, scars, and various dermatologic lesions.



Collagen Shaping lifting and firming

A unique technique that firms, fills and lifts the skin from inside.

Thin acupuncture needles that are painlessly inserted in the dermal layer, allows the doctor to make new skin supportive thread from patient's own collagen. A light, imperceptible electrical impulse from the needle produces vibration of the tissue ending with strong contraction of the skin collagen sticking to the needle. In this way the doctor can "wind up" the relaxed and sagging collagen fibers into a firm strand.

These strands form a supportive network making the skin tighter and lifted. At the same time all layers of patient's skin are stimulated, resulting in extra new young collagen, elastin and hyaluronic acid. The skin condition is greatly improved, and aging is slowed down.

**Collagen Shaping can be
performed using Needle Shaping**

LipoPlexr fat reduction

As the skin ages, unwanted localized fat deposits often develop spontaneously. Think of areas such as above and below the eyes, nose-lip fold, "hamster cheeks" and double chin. These are generally difficult and only in an invasive way to remove or relieve.

LipoPlexr delivers a lipolytic (fat-depleting) pulse locally and specifically. Small and thin acupuncture needles are inserted painlessly into the area concerned, after which a barely noticeable, slight electrical impulse is delivered to the fat cell. This creates a process of fat cell breakdown at a natural pace for the body and without damage on the surrounding tissue. Surgical fat removal or introduction of lipolytic injections can thus be avoided.

Anesthesia is not necessary for this LipoPlexr treatment and there is no down time. Depending on the size and compactness of the zone, between two and five treatments are required on average.


**Lipo Plexr can be performed
using OFF**

Endo Plexr®

Different technologies of the Omnia Skin Systems are developed to connect perfectly with each other. EndoPlexr is one of them which is one of the most powerful combination treatments frequently recommend by the doctors. This method is scientifically recognized in literature as HLH® technique of Dr. Sotiris Tsioumas, one of the most experienced Doctors in minimally invasive surgery (Soft Surgery) with GMV Medical Devices.

EndoPlexr is a hybrid method in which the effect of the Plexr and the CollagenShaping reinforce each other significantly. An EndoPlexr treatment simultaneously reduces excess and sagging skin and creates support in the skin for lifting and filling.





How does it work?

The first step is to “wind up” the collagen present in your skin into the threads. These threads of own collagen are relocated by the doctor at strategic facial positions to form a support network to fill and lift the problematic skin areas.

These threads are then fixed and strengthened by the energy impulses of the Plexr. At the same time, the doctor treats the wrinkles and excess sagging skin with Plexr further reducing and tightening it.


The EndoPlexr offers a fantastic multi solution for most of the skin needs.

In addition to the immediately visible aesthetic results, unlike the most aesthetic methods offers the EndoPlexr enormous anti-aging benefits in the long term as well.

The impulses delivered to the skin during treatment cause intense bio-stimulation. As a result, the overall skin condition and quality of your skin improves from the inside out and the skin regains its elasticity and shine. The patient will experience an intense firmness and rejuvenation due of the creation of new collagen.


Due to the harmonization of the slowed skin metabolism and the normalization of cell division, the age-diluted and sagging skin shows 4-6 weeks after the treatment an impressive tightness and volume improvement of no less than 58%. The excellent and very natural aesthetic results, combined with the high safety of the long-lasting bio-stimulation, makes the EndoPlexr so special and absolutely unique approach in the esthetic world.

Providing
a unique
multi-level action
in a **single** treatment





PLEXR OMNIA
the beauty
in your **hands**



Results




Before / After 8
months
Collagen Shaping
& Lipo Plexr




Before / After 9
months
Collagen Shaping /
Lipo Plexr




Before / After 3
months
Collagen Shaping




Before / After 1
month
Endo Plexr



Before / After 6
months
Collagen Shaping



Before / After 1
month
Endo Plexr



Credits: Anna Malecka (Poland) - Be Fit & Beautiful / Branko Mitrovic (Nederland) & Alina Markou (Cyprus) - Plexr Nederland
All rights are reserved. Reproduction is strictly forbidden.





www.gruppogmv.com

